COMPETENCY: Blood Pressure Measurement



Million Hearts Blood Pressure Project Parish Nurse:

Performance	Method of Evaluation	Met	*Unmet	Evaluator	Parish
Steps				Date &	Nurse,
				Initials	Date &
					Initials
View video clip at:					
http://www.nejm.org/doi/full/10.1056/nejmvcm08001	View 100% of video clip				
<u>57</u> .					
Describe how to determine the correct cuff size for					
the patient					
 Identify rubber bladder in cuff 	Return Demonstration				
• Bladder covers 80% circumference of arm					
• Width of bladder 40% circumference of arm					
Demonstrates correct patient positioning					
 Seated with feet flat on the floor 					
• Leaning against back of the chair not on arm	Return Demonstration				
Arm fully supported, brachial artery at heart					
level					
• Upper arm bare, do not apply cuff over clothing					
Demonstrates correct cuff placement					
 Apply cuff snugly 	Return Demonstration				
Center of bladder is directly above brachial					
artery					
Bottom edge of cuff is 1" above elbow					
Demonstrates palpatory technique of estimating BP					
 Inflate manometer while palpating radial pulse 	Return Demonstration				

•	Note level at which radial pulse disappears					
•	Release air from cuff slowly					
	Performance	Method of Evaluation	Met	*Unmet	Evaluator	Parish
	Steps				Date and	Nurse,
	-				Initials	Date and
						Initials
Der	monstrates correct technique of measuring BP					
•	Place stethoscope on brachial artery	Return Demonstration				
•	Inflate cuff 20-30 mm Hg above the point where	or				
	radial pulse disappeared	Verification with instructor				
•	Systolic reading – first two consecutive sounds	using a teaching				
•	Diastolic reading – at the level where the sound	stethoscope				
	disappears, 2 mm Hg below the last sound					
Verbalizes 2 things that may alter BP readings						
•	Caffeine • Smoking					
•	Medications • Pain					
•	Incorrect Cuff • Anxiety	Verbal Test				
•	Bad Equipment • Talking					
•	Bad Pt. Positioning • Needing to Urinate					
•	Incorrect Technique • Nurse Bias					
•	Level of Brachial Artery					
Verbalize the national standard for the normal BP reading for an adultsystolic <120 diastolic <80		Verbal Test				

* It performanc	e step unmet, w	vrite action plan be	elow:		
Evaluator:					

